|           |   |    |
|-----------|---|----|
| English   | Instructions for use: Dental Prosthetics        | 3  |
| Español   | Instrucciones de uso: Prótesis Dentales         | 6  |
| русский   | Инструкция по применению: Зубные протезы        | 9  |
| Deutsch   | Gebrauchsanweisung: Dentale Prothetik           | 12 |
| Français  | Instructions d'utilisation: Prothèses Dentaires | 15 |
| Italiano  | Istruzioni per l'uso: Protesi Dentali           | 18 |
| Português | Instruções de utilização: Próteses Dentárias    | 21 |
| Türkçe    | Kullanım talimatları: Dental Protetik Ürünler   | 24 |
| 简体中文      | 使用说明: 牙科义齿                                      | 27 |
| 日本語       | 使用説明書: 歯科プロテーゼ                                  | 30 |
| 한국어       | 적응증: 치과 보철                                      | 33 |
| عربي      | إرشادات الاستخدام: البدليات السنية              | 36 |
| Polski    | Instrukcja obsługi: Protetyka Dentystyczna      | 39 |

manufactured by:  
BioHorizons

2300 Riverchase Center  
Birmingham AL, 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
www.biohorizons.com

EU Authorised Representative:  
QUALITY FIRST INTERNATIONAL  
Suites 317/318 Burford Business Centre  
11 Burford Road, Stratford  
London E15 2ST United Kingdom  
Telephone +44-208-221-2361  
Telefax +44-208-221-1912

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

| Symbol  | Symbol Description  |
|---|---|
|    | Caution   |
|    | Electronic instructions for use   |
|    | Manufacturer  |
|    | BioHorizons products carry the CE mark and fulfill the requirements of the Medical Devices Directive 93/42/EEC                  |
|   | Reference / article number  |
|  | Lot/ batch number   |
|  | Do not re-use   |
|  | Use-by-date   |
|  | Sterile by gamma irradiation  |
|  | Date of manufacture   |
|  | Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician |
|  | Home position   |

This document supersedes all prior revisions. Original language is English.

## DESCRIPTION

BioHorizons prosthetic components are intended for the restoration of dental implants within the specific indications of each implant system. The label on each prosthetic package contains important product information including whether the prosthetic is supplied sterile or non-sterile.

## INDICATIONS FOR USE

Tapered Internal and Internal Prosthetics are intended for the restoration of Tapered Internal and Internal dental implants and Zimmer® Dental ScrewVent® and Tapered ScrewVent® dental implants. Laser-Lok 3.0 Prosthetics are intended for the restoration of BioHorizons Laser-Lok 3.0 dental implants. Single-stage Prosthetics are intended for the restoration of BioHorizons Single-stage dental implants. External Prosthetics are intended for the restoration of BioHorizons External dental implants. Overdenture Prosthetics are intended for the restoration of BioHorizons Overdenture dental implants. One-piece Prosthetics are intended for the restoration of BioHorizons One-piece dental implants. Ceramic abutments are only intended for cement-retained cases in the anterior esthetic zone. BioHorizons Angled Multi-unit Abutments contain a unique LIGHT BLUE abutment screw specifically designed for use in those abutments. Multi-unit and Abutment for Screw prosthetics are intended for use as a support for multiple-implant supported prostheses. Ball Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible (3.0mm Ball Abutments being intended for use with 3.0mm dental implants). OD Secure Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible and maxilla. Easy Ti Temp abutments are intended for long term temporary restorations (>30 days); hexed for single-unit and non-hexed for multiple-unit restorations. PEEK prosthetics are intended for short term temporary restorations (<30 days).

## CONTRAINDICATIONS

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. BioHorizons Multi-unit and Abutment for Screw prosthetics are contraindicated for single tooth replacements.

## DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions For Use that accompany BioHorizons products. Healing abutments should be hand tightened (10-15 Ncm). Abutment screws must be torqued to 30 Ncm to prevent screw loosening.

| Platform                  | Zimmer Tapered ScrewVent Implants* | Zimmer ScrewVent Implants* |
|---------------------------|------------------------------------|----------------------------|
| 3.5mm Internal Connection | TSVBx<br>TSV4Bx<br>TSVT4Bx         | TSVHx<br>TSV4Hx            |
| 4.5mm Internal Connection | TSVWBx<br>TSVTWBx<br>TSVWHx        | SVWBx<br>SVWHx             |
| 5.7mm Internal Connection | TSV6Bx<br>TSVT6Bx<br>TSV6Hx        | N/A                        |

\*Where variable x = implant length.

Secure Custom Cast Abutments to the supplied polishing protector to protect the implant/abutment interface during laboratory procedures.

Single-stage Titanium Abutments for Laser Weld should be secured to the working model by torquing the accompanying screw to 30Ncm prior to laser welding. Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling. BioHorizons Multi-unit Prosthetic Screws must be torqued to 15 Ncm to prevent screw loosening. Multi-unit Abutment screws must be torqued to 30 Ncm.

#### **WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic components. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding specific IFU.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment. Use of the LIGHT BLUE abutment screws (packaged with Angled Multi-unit Abutments) with incompatible abutments may result in abutment and/or abutment screw failure.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

#### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) abutment screw and/ or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/ or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

## **HANDLING AND STERILIZATION**

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. If the product is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined and the surgical site has been prepared. Cleaning of Laser-Lok abutments after preparation - In addition to standard abutment cleaning protocols, clinicians should ultrasonically clean the abutment with Enzymax or equivalent for a minimum of two (2) minutes prior to placement in the restorative site. Ultrasonic cleaning is required to remove particulate created during abutment preparation. If desired, a prosthetic supplied non-sterile or a prepped abutment may be sterilized - place product in an approved sterilization bag or wrap and run through one of the following qualified sterilization cycles: Non ceramic prosthetics - 1. Prevacuum Steam: 132°C (270°F) for five (5) minutes minimum. 2. Gravity Steam: 121°C (250°F) for twenty (20) minutes minimum. Ceramic prosthetics - dry heat sterilize at 170°C (338°F) for sixty (60) minutes. Do not steam sterilize ceramic abutments.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc. Z-CAD® is a registered trademark of Metoxit AG.

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