

Tapered Short Drill Usage Chart

Part Description	Part Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Starter Drill, 2.2/3.2mm	TDS32HD																				
4.6 Intermediate Drill, 2.8/3.65mm	TDS37HD																				
4.6 Final Drill, 4.1mm	TDS41HD																				
5.8 Intermediate Drill, 3.6/4.65mm	TDS47HD																				
5.8 Final Drill, 5.4mm	TDS54HD																				

TSG 4.6mm Drill, #1, 2.5mm	TSGG1																				
TSG 4.6mm Drill, #2, 2.2/3.2mm	TSGG2																				
TSG 4.6mm Drill, #3, 3.7mm	TSGG3																				
TSG 4.6mm Drill, Dense Bone, 4.1mm	TSGGDB																				
TSG 5.8mm Drill, #1, 2.5mm	TSGB1																				
TSG 5.8mm Drill, #2, 2.2/3.2mm	TSGB2																				
TSG 5.8mm Drill,, #3, 4.7mm	TSGB3																				
TSG5.8mm Drill, Dense Bone, 5.4mm	TSGBDB																				

Important Note:
All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.



customer care: 888.246.8338
shop online at store.biohorizons.com