



Tapered Internal HD surgical kit (TSK4000)

Includes the instrumentation required to place:
Tapered Plus, Tapered Tissue Level, Tapered 3.0 and
Tapered Internal including the 3.4mm diameter and 18mm lengths.



Straight Parallel Pins (2 per kit) 144-100

20° Angled Parallel Pins (2 per kit) 144-200

3.0 Implant-level Driver, Handpiece, Regular TP3IDHR

3.5/4.5mm HD Implant-level Driver, Handpiece TYGIDH

5.7mm HD Implant-level Driver, Handpiece TBIDH

Internal Abutment-level Driver, Handpiece PHA

3.0 Implant-level Driver, Ratchet, Regular TP3IDRR

3.5/4.5mm HD Implant-level Driver, Ratchet TYGIDR

5.7mm HD Implant-level Driver, Ratchet TBIDR

Internal Abutment-level Driver, Ratchet PRA

Ratchet 130-000

Spare slots can be used to hold preferred surgeon instrumentation such as stop drills and extended shank drills. The implant staging area is designed to temporarily secure implant vials during surgery only. Do not attempt to sterilize implant vials in the surgical tray.

customer care: 888.246.8338 • shop online at store.biohorizons.com

ARCHIVED DOCUMENT

BIOHORIZONS®

Tapered Internal HD
surgical kit (TSK4000)



L02025



REV C MAY 2017